





# Safety Lancet Cross-Reference Guide

## Competitors Lancets





## SurgiLance™/SurgiLance™ Lite

BRAND/ MANUFACTURER	PRODUCT	CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH	PRODUCT/ CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH
ACCU-CHEK SAFE-T-PRO (ROCHE)  	Safe-T-Pro	951	Blue/ White	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
ACCU-CHEK SAFE-T-PRO PLUS (ROCHE)  	Safe-T-Pro Plus	BID448622001	Purple/ White	Needle	23	1.3	SLN100	Yellow	Needle	21	1.0
						1.8	SLN200	Gray	Needle	21	1.8
						2.3	SLL21G180	Purple	Needle	21	1.8
							SLN240	Orange	Needle	21	2.2
							SLL21G220	Lime Green	Needle	21	2.2
ACTI-LANCE (HTL STREFA)  	Lite	7141	Purple	Needle	28	1.5	SLL28G140	Yellow	Needle	28	1.4
	Universal	7142	Blue	Needle	23	1.8	SLN200	Gray	Needle	21	1.8
	Special	7143	Gold	Needle	17	2.0	SLL21G180	Purple	Needle	21	1.8
							SLB200	Green	Blade	18	1.8
BD MICROTAINER  	Low Flow	366592	Purple	Needle	30	1.5	SLL28G140	Yellow	Needle	28	1.4
	Medium Flow	366593	Pink	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
	High Flow	366594	Blue	Blade	–	2.0	SLL21G180	Purple	Needle	21	1.8
							SLB200	Green	Blade	18	1.8

# Safety Lancet Cross-Reference Guide

## Competitors Lancets





## SurgiLance™/SurgiLance™ Lite

BRAND/ MANUFACTURER	PRODUCT	CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH	PRODUCT/ CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH
<b>CARDINAL HEALTH</b>  	Low Flow	SLLF100/200	Lavender	Needle	28	1.8	SLL28G180	Gray	Needle	28	1.8
	Medium Flow	SLMF100/200	Yellow	Needle	23	1.8	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
	High Flow	SLHF100/200	Orange	Needle	21	2.0	SLN240	Orange	Needle	21	2.2
	Dual Flow	SLDF100/200	Light Blue	Needle	18	1.8	SLL21G220	Lime Green	Needle	21	2.2
						SLB200	Green	Blade	18	1.8	
<b>ERGOLANCE (HTL STREFA)</b>  	Micro Flow	8484/8491	Green	Needle	30	–	SLL28G140	Yellow	Needle	28	1.4
	Normal Flow	8495/8492	Blue	Needle	25	–	SLL28G180	Gray	Needle	28	1.8
	High Flow	8486/8493	Pink	Needle	21	–	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
<b>EZ-LANCE (AMBIMEDINC)</b>  	Lite	ELY100	Yellow	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
	Regular	ELG100	Green	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
	Regular Plus	ELO100	Orange	Needle	21	2.2	SLN240	Orange	Needle	21	2.2
						SLL21G220	Lime Green	Needle	21	2.2	
Ultra	ELR100	Blue	Blade	21	2.2	SLB200	Green	Blade	18	1.8	
<b>GENIE (BD)</b>  	Low (Single Drop)	366579	Purple	Needle	28	1.25	SLL28G140	Yellow	Needle	28	1.4
							SLN100	Yellow	Needle	28	1.0
	Low (Single Drop)	366583	Yellow	Needle	23	2.25	SLN240	Orange	Needle	21	2.2
	Low	366580	Pink	Blade	–	1.0	SLL21G220	Lime Green	Needle	21	2.2
	Medium	366581	Green	Blade	–	1.5	SLB200	Green	Needle	18	1.8
	Medium-High	366582	Blue	Blade	–	2.0	SLB250	Blue	Blade	18	2.3

# Safety Lancet Cross-Reference Guide

## Competitors Lancets




## SurgiLance™/SurgiLance™ Lite

BRAND/ MANUFACTURER	PRODUCT	CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH	PRODUCT/ CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH
<b>HAEMOLANCE PLUS</b> (HTL STREFA)  	Micro Flow	7370/7360	Light Blue	Needle	28	1.6	SLL28G140	Yellow	Needle	28	1.4
	Low Flow	7371/7361	Blue	Needle	25	1.4	SLL28G180	Gray	Needle	28	1.8
	Normal Flow	7372/7362	Green	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
	High Flow	7373/7363	Yellow	Needle	18	1.8	SLN300	Pink	Needle	21	2.8
							SLL21G280	Pink	Needle	21	2.8
	Max	7374/7364	Purple	Blade	–	1.6	SLB200	Green	Blade	18	1.8
<b>MEDISAFE SOLO</b> (HTL STREFA)  	–	7130	Green	Needle	29	1.5	SLL28G140	Yellow	Needle	28	1.4
	–	7140	Pink	Needle	23	2.0	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
<b>MEDLANCE PLUS</b> (HTL STREFA)  	SuperLite	7241/7242	Lime	Needle	30	1.2	SLL28G140	Yellow	Needle	28	1.4
	Lite	7043/7243	Lavender	Needle	25	1.5	SLL28G180	Gray	Needle	28	1.8
	Universal	7044/7244	Blue	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
	Extra	7045/7245	Green	Needle	21	2.4	SLN240	Orange	Needle	21	2.2
						SLL21G220	Lime Green	Needle	21	2.2	
	Special	7046/7246	Yellow	Blade	–	2.0	SLB200	Green	Blade	18	1.8
<b>MEDLINE</b> PUSH BUTTON  	–	MPHSAFETY28	Lavender	Needle	28	1.8	SLL28G180	Gray	Needle	28	1.8
	–	MPHSAFETY21	Orange	Needle	21	2.2	SLN240	Orange	Needle	21	2.2
							SLL21G220	Lime Green	Needle	21	2.2

# Safety Lancet Cross-Reference Guide

## Competitors Lancets





## SurgiLance™/SurgiLance™ Lite

BRAND/ MANUFACTURER	PRODUCT	CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH	PRODUCT/ CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH
<b>MEDLINE PRESSURE ACTIVATED</b>  	—	MPHST28	Light Blue	Needle	28	1.6	SLL28G140	Yellow	Needle	28	1.4
	—	MPHST21	Green	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
<b>MEDLINE SIDE ACTIVATED</b>  	—	MPHSTSIDE28	Purple	Needle	28	1.8	SLL28G180	Gray	Needle	28	1.8
	—	MPHSTSIDE21	Pink	Needle	21	2.4	SLN240	Orange	Needle	21	2.2
							SLL21G220	Orange	Needle	21	2.2
<b>NOVLANCE (CH DIAGNOSTICS)</b>	Novlance	4225	Blue	Needle	25	1.8	SLL28G180	Gray	Needle	28	1.8
<b>PROLANCE (HTL STREFA)</b>  	Micro Flow	7599	Light Blue	Needle	28	1.6	SLL28G140	Yellow	Needle	28	1.4
	Low Flow	7594	Blue	Needle	25	1.4	SLL28G180	Gray	Needle	28	1.8
	Normal Flow	7601	Green	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
								SLL21G180	Purple	Needle	21
	High Flow	7620	Yellow	Needle	18	1.8	SLN300	Pink	Needle	21	2.8
	Max Flow	7609	Purple	Blade	—	1.6	SLL21G280 SLB200	Pink Green	Needle Blade	21 18	2.8 1.8

# Safety Lancet Cross-Reference Guide

## Competitors Lancets

## SurgiLance™/SurgiLance™ Lite

BRAND/ MANUFACTURER	PRODUCT	CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH	PRODUCT/ CATALOG #	COLOR	NEEDLE/ BLADE	GAUGE	DEPTH
<b>SAFE-T-LANCE</b> (SMITHS MEDICAL)  	Micro Flow	1025	Light Blue	Needle	28	1.8	SLL28G180	Gray	Needle	28	1.8
	Low Flow	1020	Blue	Needle	25	1.8	SLLG28220	Orange	Needle	28	2.2
	Normal Flow	1014	Green	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
	High Flow	1008	Yellow	Needle	18	1.8	SLL21G180	Purple	Needle	21	1.8
							SLN300	Pink	Needle	21	2.8
	Maximum Flow	1005	Purple	Blade	–	–	SLL21G280	Pink	Needle	21	2.8
						SLB250	Blue	Blade	21	2.3	
<b>STAT AUTO</b> (STAT MEDICAL)  	Lite	SAL-100	Purple	Needle	28	1.6	SLL28G140	Yellow	Needle	28	1.4
	Tender	SAP-100	Blue	Needle	26	1.8	SLL28G180	Gray	Needle	28	1.8
	Comfort	SAT-100	Pink	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
	Normal	SAC-100	Gray	Needle	21	2.4	SLL21G180	Purple	Needle	21	1.8
						SLN240	Orange	Needle	21	2.2	
						SLL21G220	Lime Green	Needle	21	2.2	
<b>UNISTIK 3</b> (OWEN MUMFORD)  	Comfort	AT104	Lavender	Needle	28	1.8	SLL28G180	Gray	Needle	28	1.8
	Normal	AT100	Yellow	Needle	23	1.8	SLN200	Gray	Needle	21	1.8
	Extra	AT101	Orange	Needle	21	2.0	SLL21G180	Purple	Needle	21	1.8
							SLN240	Orange	Needle	21	2.2
						SLL21G220	Lime Green	Needle	21	2.2	
<b>VITALCARE AUTO</b>  	Junior	FTM065	Pink	Needle	26	1.8	SLL28G180	Gray	Needle	28	1.8
	Soft	FTM066	Blue	Needle	21	1.8	SLN200	Gray	Needle	21	1.8
							SLL21G180	Purple	Needle	21	1.8
	Normal	FTM067	Beige	Needle	21	2.4	SLN240	Orange	Needle	21	2.2
	Extra	FTM069	Green	Needle	21	3.0	SLL21G220	Lime Green	Needle	21	2.2
						SLL21G280	Pink	Needle	21	2.8	